





HOW TO RECOGNIZE AN ORIGINAL EVOLIS HIGH TRUST® CONSUMABLE

Hereafter is a description of the several points to check depending on your ribbons.

RIBBONS PLACED IN CASSETTE for Primacy/Zenius printers:

• Check the Evolis High Trust® and Evolis logos on the label.



• We placed two 'O' from the Evolis logo on each plastic ribbon core.



RIBBONS W/O CASSETTE EASY4PRO N5F008xxx, N5H004NAA, N6F003xxx, N7H006NAA, NCT023NAA, N2F010NAA (for Primacy printers):



RIBBONS WITH PART NUMBERS R3012, R3011, R3511, R3411, R3013, R3314, R3514, IN INDIVIDUAL BOXES (for Tattoo/Pebble/Dualys/Securion/Quantum printers):

- Check the Evolis High Trust® and Evolis logos on the box.
- We placed a **holographic sticker** on the ribbon box, with several security features, such as a rainbow color fine lines guilloche design, nano text, two channel images...
- We also integrated a unique feature: the floating image. To check the authenticity of your ribbon, you just have to pass the sticker under a standard light and the 'O' from the Evolis logo will appear.





• On each side of the ribbon cores, check the logo



MONOCHROME RIBBONS WITH PART NUMBERS R2011, R2211, R2131, R2121 IN SEALED PLASTIC BAG (for Tattoo/Pebble/Dualys/Securion/Quantum printers):

- Check the Evolis High Trust® and the Evolis logos on the label.
- On each side of the ribbon cores, check the logo





RIBBONS for Avansia printers:

 \bullet Check the Evolis High $\textbf{Trust}^{\textcircled{\tiny{0}}}$ and Evolis logos on the label.





• Check the **Evolis logo** on the RFID tag.

